

**IMMUNOSTICS, INC.**

3505 SUNSET AVENUE  
OCEAN, NEW JERSEY 07712  
732-918-0770 800-722-7505 FAX 732-918-0618

**MATERIAL SAFETY DATA SHEET**

MSDS# IDM001  
Revised: 01/2008

NFPA/HMIS Hazard Rating Health: 1\*  
Fire: 0  
Reactivity: 0

\*Precautions required: See Sections II and IV

**SECTION I: IDENTIFICATION OF PRODUCT**

PRODUCT NAME

CATALOGUE NO.

Immuno Detector - MONO

IDM-25

**PRODUCT DESCRIPTION**

A qualitative immunoassay for the detection of infectious mononucleosis antibodies in human whole blood, serum or plasma.

**SECTION II: HAZARDOUS INGREDIENTS****REAGENTS CONTAINED IN KIT****COMPOSITION**

Diagnostic test strip in a plastic cassette  
Packaged in a foil pouch.

Combination of mouse monoclonal  
Antibodies and polyclonal antibodies.

Positive(+) and Negative(-) controls

Human serum or defibrinated plasma with  
Plasma with 0.1% Azide as preservative.

Stabilized Buffer.

PBS pH7.2±0.2 with 0.1% Azide as  
preservative.

All components derived from human source materials have been tested and found to be non-reactive for HbsAg and antibodies to HIV and HCV. Given that no known test offers complete assurance that infectious agents are absent, all materials derived from human blood should be handled as if capable of transmitting infection.

Post-It® Fax Note	7671	Date	9/23	# of pages	3
To	Tarheel	From	Cindy		
Co./Dept.	Carmen	Co.	Immunostics		
Phone #		Phone #	732-918-0770		
Fax #	910-762-6022	Fax #			

The information provided in this Material Data Safety Sheet has been compiled from our experience and data presented in various technical publications. It is believed to be correct, however, it is the user's responsibility to determine the suitability of this information for the adoption of safety precautions deemed necessary. Immunostics, Inc. shall not be held liable for any damage resulting from handling or from contact with the product described in this MSDS. We reserve the right to update MSDS sheets from time to time as new information becomes available. It is the responsibility of the user to verify that they have the latest revision available.

**SECTION III: PHYSICAL AND CHEMICAL CHARACTERISTICS**

Chemical Family: Blood serum protein

Fire and Explosion Hazards: None

Reactivity: Sodium azide preservative may react with lead and copper plumbing to form highly explosive metal azides. When disposing, flush with a large quantity of water to prevent azide buildup.

Incompatibility: No known materials that must be avoided.

Toxicity and Carcinogenicity: No OSHA exposure limit has been established and no carcinogenic effect has been reported. No known toxic material at a concentration of 1% or higher.

**SECTION IV: HEALTH HAZARDS AND HANDLING PRECAUTIONS****PRIMARY ROUTES OF ENTRY****PRECAUTIONS**

**Ingestion:** Avoid hand-to-mouth contact when handling human source materials. Wash hands thoroughly after handling, even when gloves have been worn. Do not eat, drink or apply cosmetics in the area where human source material is handled. Do not pipet by mouth.

**Skin:** Wear gloves and especially cover any cuts, abrasions or skin lesions. Dispose of gloves pipettes, stirrers, test cards and used reagent containers as biohazardous material. Wash hands thoroughly after removing gloves. Use extreme caution with any sharp object to avoid percutaneous exposure to human source material. Wear outer protective garments such as a lab coat or gown.

**Inhalation:** If splash or aerosol can be created from manipulation of this product, use a surgical mask or similar respiratory protection to cover nose, mouth and mucous membranes.

**Eyes:** If splash or aerosol can be created from manipulation of this product, use chemical safety goggles, face shield or splash shield as appropriate to prevent eye exposure.

The information provided in this Material Data Safety Sheet has been compiled from our experience and data presented in various technical publications. It is believed to be correct, however, it is the user's responsibility to determine the suitability of this information for the adoption of safety precautions deemed necessary. Immunogistics, Inc. shall not be held liable for any damage resulting from handling or from contact with the product described in this MSDS. We reserve the right to update MSDS sheets from time to time as new information becomes available. It is the responsibility of the user to verify that they have the latest revision available.

**SECTION V: EMERGENCY AND FIRST AID PROCEDURES**

- Ingestion: If conscious, wash out mouth with water. Call a physician.
- Skin Contact: Wash thoroughly with soap and water. Remove contaminated clothing. Call a physician.
- Inhalation: Remove from source to fresh air. If breathing becomes difficult, call a physician.
- Eye Contact: Flush with large amounts of water or sterile eye wash. Use fingers to separate the eyelids for effective flushing. Call a physician.

**SECTION VI: SPILL/LEAK PROCEDURES**

All materials derived from human blood should be handled as if capable of transmitting infection. Other materials that are contaminated with human source material should be considered to be capable of transmitting infection.

If material is released or spilled, wear all appropriate protective equipment described in Section IV before cleaning up the spill or handling contaminated material. Wipe up the spill and dispose of the contaminated materials in a bag to be identified and treated as biohazardous waste. Avoid creating aerosols or dust while cleaning up a spill.

After the material has been picked up and contained in a bag, wash the spill site with a disinfectant cleaner.